

geemarc DALLAS 10

User guide

T.B.A.

If you are using this telephone with a hearing aid, set your hearing aid to the «T» setting

T e l e p h o n e

SETTING UP

T.B.A.
Drawing description



(*) Classified TNV-3 according to EN60950 standard.

1 - CONNECTING TELEPHONE

- 1 - Connect the handset cord.
- 2 - Connect the line cord to the socket located at rear of telephone then connect the PTT plug*.

SETTING UP

2 - RINGER SETTINGS

A switch located on right side allows the ringer to be switched On () or Off (.

3 - WALL MOUNTING

To wall-mount your telephone, slide out the clip turn around and slide back into the wall mounting position (located in the base unit, just below the ear piece). This will ensure that the handset stays firmly in position. Make 2 holes in the wall 85mm apart, insert the wall plugs and the wall mounting screws. The screws should protrude from wall by 6-7 mm. Place the phone onto the screw-heads and slide down to secure.

4 - TONE/PULSE SWITCH (T/P)

The Tone/Pulse switch is factory pre-set to Tone position (**T**). In the UK, all telephone exchanges now use Tone dialing. If your phone does not dial out it is probably being used from an older private switch board (PBX). In this case, the Tone/Pulse switch must be moved to Pulse (**P**).

USING THE PHONE


1 - ANSWERING INCOMING CALL

- When an incoming call is received the phone will ring and flash. To answer the call lift the handset and speak.
- On completion of the call, carefully replace the handset in the cradle to release the line.

2 - VOLUME RECEIVING ADJUSTMENT

You can increase the volume according to your own hearing requirements.

Use the **VOL** button, located on the right side, to adjust the level according to your own hearing requirements.





USING THE PHONE

3 - MAKING A CALL


- Lift the handset.
- Await the dial tone and dial required number.
- On completion of the call, carefully replace the handset in the cradle.

4 - SENDING LEVEL ADJUSTMENT

If your voice is quiet or called parties have difficulty hearing you, you can increase the sending level volume with the switch  located on the right side.



5 - LAST NUMBER REDIAL

- Lift the handset.
- Await the dial tone and press  button.
- The previously dialled number will be automatically redialled (not in use for memory numbers).



USING THE PHONE

6 - R, * & # BUTTONS

These are used with the new services provided on digital exchanges. For details please contact your network operator.

7 - STORING TELEPHONE NUMBERS

Your phone can store 3 direct and 10 indirect memories (16 digits maximum for each memory).

- Lift the handset.
- Press  button.
- Dial the telephone number you wish to store.
- Press  button.
- Press one of the direct memory buttons, **M1**, **M2**, **M3**.

Or press **MEM** and any button **0** to **9** on the key pad under which you wish to store the telephone number in indirect memory.

- Replace the handset.

USING THE PHONE

Note : Numbers will be lost if you disconnect your phone. Entering a new number in memory will automatically erase the previous number.

8 - DIALLING STORED NUMBER

- Lift the handset and await the dial tone.
- Press appropriate direct memory button **M1**, **M2**, or **M3**. The stored number will be dialled out automatically.

Or press the **MEM** button followed by the keypad **0** to **9** button under which the desired number is stored. The stored number will be dialled out automatically.

To register your stored phones numbers, use the index located on the base.

GUARANTEE

From the moment your Geemarc product is purchased, Geemarc guarantee it for the period of one year. During this time, all repairs or replacements (at our option) are free of charge. Should you experience a problem then contact our helpline or visit our website at www.geemarc.com. The guarantee does not cover accidents, negligence or breakages to any parts. The product must not be tampered with or taken apart by anyone who is not an authorised Geemarc representative. The Geemarc guarantee in no way limits your legal rights.

IMPORTANT: YOUR RECEIPT IS PART OF YOUR GUARANTEE AND MUST BE RETAINED AND PRODUCED IN THE EVENT OF A WARRANTY CLAIM

Please note: The guarantee applies to the United Kingdom only.

Declaration : Hereby Geemarc Telecom SA declares that this product is in compliance with the essential requirements and other relevant provisions of the Radio and Telecommunications Terminal Equipment Directive 1999/5/EEC and in particular article 3 section 1a, 1b and section 3.

Telephone connection : Voltages present on the telecommunication network are classified TNV-3 (Telecommunication Network Voltage) according to the EN60950 standard.



**For product support and help visit our website
at www.geemarc.com
telephone 01707 384438
or fax 01707 372529**